

# SELAND CARE

THERMAL SUIT  
FOR PEOPLE  
WITH DISABILITES



It is a thermal neoprene suit, specially designed to be easily put on and taken off and that adapts to each individual. It maintains the body temperature, protects from shocks and regulates the balance of a person with reduced mobility or people with moderate, severe or total disability. It allows you to enjoy the hydrotherapy session with only one companion or professional.

The time it takes to put it on can be between three and six minutes, depending on the degree of disability each person has.

**ALWAYS SUPERVISED BY A COMPANION WHO CAN ASSIST IN CASE OF NEED (avoid turning around).**

# THERMAL SUIT FOR PEOPLE WITH DISABILITIES

Thermal protection to stay four times longer in the water, abrasion protection and the blows, with a system that regulates the balance with negative or positive buoyancy according to the needs.

People with disabilities are the biggest candidates to use this kind of suits:

- When they are in the water their body is freed from these physical bonds, which leads to paralysis and enjoy weightlessness.
- Water is one of the most suitable means of rehabilitation and physical maintenance for people with disabilities. The use of a thermal suit is very important, since they do not move enough in the water to generate heat and the temperature drops quickly. With the neoprene we manage to keep the body heat longer and they can do the rehabilitation exercises without having to stop them immediately, due to the shivering produced by the cold.
- The neoprene has great protection against all impacts and rubbing especially for people with disabilities, who do not control their body movements well, so they are much more vulnerable to injury.
- The body is balanced by the positive buoyancy of the neoprene and the negative buoyancy of a few individual weights, so that it helps us to keep the person in the proper posture in the water, for their rehabilitation.
- In addition, he wears a floating collar that keeps his head out of the water for those who need it, which can be adjusted by adding or changing the position of the inner reinforcements.

All these characteristics make these people free from their body, floating in the water, enjoying the moment and helping both caregivers and rehabilitation professionals by facilitating their work in the pool.

In this way, these people can be easily maintained with the support of one person.



# SUIT PLACEMENT

It is considered an easy to wear system, since the suit opens up whole, leaving only one piece. In this way the suit is placed on a flat surface (stretcher or table).

The velcro is opened leaving the inside, which is the soft part in sight, so we can place the person lying on top.  
(Photo 1)



It's closing and wrapping the person:  
Start at the chest, attaching the Velcro on the two front legs.  
(Photo 2)



Then, close their legs, starting at the ankles and going up the thighs, to the groin.  
(Photo 3)




Then it is followed by the arms, starting from the wrists to the neck.  
(Photo 4)




The collar is attached to the suit with Velcro and loops around the neck of the suit. It can also be adjusted by changing or adding foams where needed.

**CLOSE THE VELCRO AFTER YOU TAKE OFF THE SUIT.  
Velcro can damage the inside of the suit and the seams.**

# ADJUSTMENT IN THE WATER



POSITIVE FORCE  
NEOPRENE - FOAM



NEGATIVE FORCE  
NEOPRENE - FOAM

ON THE ANKLES:  
Replace the foam with the weights  
(remove all foam from thighs)



NEGATIVE FORCE  
OR WEIGHT



SELAND NEOPREN S. L.  
Antsonekoa kalea, 2 (Pol. Igeltzera. Naves 4-5)  
48610 Urduliz (Bizkaia)  
Tlf. (+34) 94 676 48 75 • [selandcare@seland.com](mailto:selandcare@seland.com)